

ACCELERATING MICROBIAL DETECTION IN CONSUMER GOODS USING RAPID STERILITY TESTING

INTRODUCTION

Ensuring sterility in consumer goods is essential to maintaining product safety and quality. However, traditional microbial detection methods, such as culture-based assays, suffer from significant limitations, particularly in time-sensitivity and detection thresholds. These standard techniques often require several days to yield results, during which production timelines are delayed, costs rise, and the risk of distributing contaminated products increases. For industries embracing just-in-time production models, such delays are not just costly, but unsustainable¹.

Conventional sterility testing techniques also struggle to detect low contamination levels, particularly in complex matrices where product ingredients and preservatives interfere with microbial detection. These factors further complicate the process, often leading to false negatives or delayed identifications. Given the increasing demand for rapid turnaround in manufacturing, there is a critical need for sterility testing solutions that offer high sensitivity without the protracted waiting times of conventional methods².

Microcalorimetry has emerged as a powerful solution for rapid sterility testing, enabling the detection of microbial contaminants with unprecedented speed and sensitivity. Unlike traditional culture-based methods, which rely on observable colony growth, microcalorimetry measures the heat produced by the metabolic activities of microorganisms. This technique provides accurate detection within hours, even in highly complex sample matrices, with a limit of detection down

to single digits CFU. The rapidity and precision of microcalorimetry have already proven effective in the pharmaceutical sector, specifically for the detection of bacterial contamination in cell and gene therapy products. The same benefits are now being explored for consumer goods, where fast, reliable, and non-destructive testing is critical for reducing product spoilage and optimizing supply chain efficiency³.

By significantly reducing the time required to confirm sterility and increasing sensitivity in low bacterial load environments, microcalorimetry provides a robust, reliable alternative to traditional testing methods. It is poised to revolutionize quality control processes, particularly in industries where swift microbial detection is essential for product safety, regulatory compliance, and minimizing production waste.

PROCEDURE

The product under investigation was mixed 1:10 in TSB media (W/V), after which 500 μ L was transferred to a vial and sealed. Thirty-two (32) samples at a time were evaluated for microbial activity in real-time in the calScreener at 37°C. A rise in metabolic activity of 5 μ W was deemed positive, and the time-to-activity was noted. Here, we tested three big-brand consumer products: liquid hand wash, shower cream, and toothpaste. To test for time-to-activity of *Escherichia coli* in the presence of the product samples were diluted and spiked into the product to a final concentration of 10 CFU per sample. Samples were transferred to calVials and incubated in the calScreener at 37°C. Samples were deemed positive when 5 μ W of metabolic heat could be detected.

1. SAMPLE PREPARATION




Various consumer products mixed with TSB 1:10

2. SAMPLE LOADING



Direct transfer to vials and spiked with ~10 CFU of *E.coli*

3. CONTINUOUS MEASUREMENT



Real-time Immediate contamination detection and time-to-activity measurements

4. REPORT



Figure 1 Testing procedure

RESULTS




Figure 2 Thermogram of products diluted in TSB with ~10 CFU of *E. coli*. The bacterial metabolism can be registered after a couple of hours in the calScreener, and after 5.5 – 6 hours the metabolic signal crosses the threshold.

In both the liquid hand wash, the shower cream, and the toothpaste, the calScreener was able to detect microbial activity within 5.5-6 hours (Fig. 2). After spiking samples with ~10 CFU of *E. coli*, we detected microbial activity less than 6 hours after initial inoculation whereas with pure TSB the detection time was ~5 hours, indicating a limited inhibitory effect of the products themselves (Fig. 3).




Figure 3 The mean time to detection of metabolic activity by *E. coli* in products diluted with TSB (n=3). The gray mark indicates the time to activity for *E. coli* in pure TSB.

CONCLUSION

Microcalorimetry offers a highly effective solution for rapid sterility testing in the consumer goods sector. Its ability to detect microbial contamination through metabolic heat production enables fast, accurate results, even in complex product matrices that often interfere with traditional methods. With a detection sensitivity as low as single-digit CFU and the capability to provide results within hours rather than days, microcalorimetry significantly enhances quality control processes. This not only reduces the risk of product spoilage but also aligns with the growing need for efficient, just-in-time production in modern supply chains. By integrating microcalorimetry, consumer goods manufacturers can ensure both product safety and operational efficiency.

REFERENCES

1. Ricardo Franco-Duarte *et al.* Advances in Chemical and Biological Methods to Identify Microorganisms—From Past to Present, Microorganisms 2019
2. Detection Limits Working Group (DLWG), Presidential Task Force on Best Practices for Microbiological Methodology, FDA
3. Christian Fricke *et al.*, Rapid Calorimetric Detection of Bacterial Contamination: Influence of the Cultivation Technique, Front. Microbiol. 2019

SYMCEL[®]